

**Product** Wireless Patient Handset

**Product Code:** NEC-05

**Description:** Wireless patient handset with Call and Reset buttons. Works with nurse call, direct to pager, and indication light

**Features:**

- Large rubber buttons for easy use by elderly patients
- Reassurance LED when transmission is sent
- Red CALL button to trigger alarm for caregiver
- Green CANCEL button will reset alarms
- Lanyard clip allows patient to connect with clothes
- Wall mount bracket for when handset is not in use
- Low battery indication

**Works with:**

- DON-20 Receiver Dongle
- ZP-01(S) Zone Protektor
- NGM-21 CordFree Monitor
- CMEX-11/21 Cordless Expander
- ILB-21 Indication Light
- PAG-11/21/31 Caregiver Alert Pager



**Technical Specifications**

Radio frequency: 433.92MHz  
Operating range: 100m (line of sight)  
Power supply: 2 x AAA alkaline batteries  
Dimensions: 115 x 42 x 25mm  
Housing material: ABS  
Net weight: 85g (without batteries)  
Operating temp.: 5°C to 40°C  
Rel. humidity: <90%

**Standards**

EMC: EN 55022  
Radio: EN 300220-1&2,  
EN 301489-1&3  
Social Alarm: EN 50130-4  
Safety: EN 60950-1  
Manufacture: ISO 9001:2008,  
ISO 13485:2003  
FCC Approved WNG-NEC-05