



Crest PermaLight™ Momentary Call Cord for Jeron® Provider® 790



Description

The PermaLight Call Cord for the Jeron Provider 790 system has LEDs that “perma”ntly light up for contiunal visibility in low-light environments.

Features

- Bright red LEDs light up pendant
- Requires no batteries
- FDA Listed

Technical Specifications

Electrical Rating:	4.6 V	Weight:	5.8 oz
Electrical Connection:	Various connection types available	Pendant Dimension:	2.75"H x 1.625"W
Button Life:	>1,000,000 cycles	Cord Length:	7'
Switch Type:	Normally Open SPST	Environment:	Operating temperature: 50° - 100°F
Warranty:	2 Years		Relative humidity: 30% to 80% (non-condensing)
Material:	UL 94V-0 Flame Rated, Latex Free PVC		Oxygen safe: Meets UL 1069 oxygen enriched atmosphere testing
Standards:	FDA Listed RoHS Compliant		Liquid ingress protection: Meets IEC 60529 IPx4 requirements
Color:	White pendant with red call button		

Crest Technical Support: 1-800-328-8908

Specifications subject to change without notice.
Visit CrestHealthcare.com for additional technical documents.