

# GripStep®

## Black XL Shoe Cover

- Non-skid in all conditions
- Safe for medical facilities
- Single use
- Non-linting
- Latex-Free



Specifications	
Construction materials	Polypropylene fabric with Anti-Skid Print
Material weights:	40 GPM (grams per square meter)
Size	Extra Large—up to size 14 Men's boots
Color	Black
Seams	Ultrasonically welded
<b>Item Number</b>	<b>Packaging</b>
GRIP4301	150 Pairs/Case



For more information visit:  
[www.baymedicalcompany.com](http://www.baymedicalcompany.com)  
 1-800-977-7888